Peritone

Instruction Manual – GB

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Manual de Instrucciones – ES

inner strength, outer confidence

Language Contents

- 4 - 33** Instructions for use - GB
- 34 - 65** Conseils d'utilisation - FR
- 66 - 96** Instrucciones de uso - ES

This instruction manual is valid for the Neen PERITONE unit. It is published by Patterson Medical. Patterson Medical does not guarantee its contents and reserves the right to improve and amend it at any time without prior notice.

Amendments will however be published in a new edition of this manual.

Revised Issue Date: October 2010 Document Number: VM-NEEN102-OM001-0

Contents

4 - 5	Warnings
6	Introduction
7	Contraindications & Precautions
8 - 9	Description of PERITONE Unit & Functions
10 - 11	Quick Start Instructions
12 - 13	Advanced Instructions
14	Recording Information
15 - 16	Statistics
17 - 18	Electrode Types & Tips
19	Care & Maintenance
20	Specifications
21	Interference
22	Troubleshooting
23 - 25	Software
26	Disposal & Applications
27 - 30	Information regarding Electromagnetic compatibility and interference (EMC)
31	Warranty
32	Clinical References

Warnings

- This unit must be used with the guidance of a Physiotherapist, Doctor or Continence Advisor.
- Type BF equipment, Continuous Operation.
- Do not insert lead wires into a mains power supply.
- Do not immerse the unit into water or any other substance.
- Do not use the PERITONE in the presence of a flammable anaesthetic gas mixture and air or with oxygen or nitrous oxide.
- If using rechargeable 9 Volt PP3 Nickel Metal Hydride batteries, be sure to use a CE approved battery charger. Never connect the PERITONE directly to a battery charger or to any other mains powered equipment.
We advise not to use Ni-Cad rechargeable batteries.
- Patient Electrodes including all skin surface electrodes, vaginal and anal electrodes are for single patient use only.
- Use only CE approved vaginal or anal electrodes i.e. PERIFORM®+ or ANUFORM®.
- Keep out of reach of children.

- To avoid the effects of electromagnetic interference never use the PERITONE device, within 4 meters of a cellular telephone or near any other powerful radio interference producing equipment that produces electrical sparks etc. In the electromyography (EMG) mode, the PERITONE device may be susceptible to strong interfering radio type emissions that may lead to temporarily increased EMG microvolt readings. The reading will immediately return to the correct value when the interference ceases. (Remember that a relaxed muscle should read below 3.5 μ volts)
- No modification of this equipment is allowed!

Symbols on the rear cabinet of PERITONE explained:

TYPE BF

Type BF
Equipment

Do not dispose in normal
dustbin (see page 26 for the
disposal instructions)

Follow
instructions
for use

Introduction

EMG (Electromyography) biofeedback is a method by which muscle activity is recorded and displayed on a unit to assist in the retraining of muscles. The PERITONE is a simple to operate, single channel EMG unit. It has been designed for use in a variety of applications, in particular for **diagnosing continence problems and for biofeedback of pelvic floor muscles**. It can also be used very effectively for Sports and General Physiotherapy applications.

The PERITONE is an accurate and sensitive single channel EMG unit that measures muscle activity down to as low as to 0.2 μV (microvolts) and up to 2000 μV . The diverse range enables the unit to measure even very weak muscle activity. It displays the biofeedback by way of bright LED lights on the front of the unit and also gives audio feedback.

As well as measuring muscle activity, it is important to follow an individual patient's progress. By locking the unit, it records the average work/rest, onset and muscle release parameters and time. All session data will be added together to form the overall average reading for the given period. This includes the number of trials performed and average muscles activity, ideal for when patients take PERITONE home and return to clinic.

Conditions known to respond to EMG:

- Incontinence • Pelvic pain • Patella femoral pain
- Neck and shoulder girdle pain • Lower Back pain
- Spinal conditions • Chronic tension headache
- Phantom limb pain • Asthma • Post-traumatic stress disorders
- Neuro feedback • Repetitive strain injury

Contraindications & Precautions


Read this operating manual before using the PERITONE

EMG should not be used:


- During menstrual period (when using a vaginal probe).
- If symptoms of bladder infection are present.
- During pregnancy (unless medically advised).
- With patients who have diminished mental capacity or physical competence who cannot handle the unit properly.
- Do not place electrodes inside mouth.
- The patient should use the unit only as prescribed.
- Do not immerse the unit in water or any other liquid.
- If in doubt about the use of the PERITONE unit, call your Doctor, Physiotherapist, Continence Advisor or your distributor for advice.

Description of PERITONE Unit & Functions

Front


Top

Rear

The concealed **LOCK** button locks the program and allows viewing of recorded information.

Quick Start Instructions

1. Insert a 9 volt PP3 Alkaline battery into the battery compartment. Can also use a rechargeable Nickel Hydride battery which has a longer life than the Ni-Cad rechargeable batteries.
2. Insert the EMG reference and EMG lead wires into the labelled sockets at the top of the unit as shown on page 9.
3. *For pelvic floor muscle biofeedback* - remove the surface electrodes and/or PERIFORM®+ or ANUFORM® from their outer packaging. For the time being, do not remove the surface electrodes from the clear plastic film.
4. Connect one surface electrode to the pin connector at the end of the EMG reference lead wire. The reference lead wire only requires one surface electrode. Attach the PERIFORM®+ or ANUFORM® cables to the pin connectors at the end of the EMG lead wire.
5. *For general muscle biofeedback* – instead of attaching the PERIFORM®+ or ANUFORM® electrode to the EMG cables attach 2 surface electrodes.
6. Make sure your skin around the treatment area is free from grease or oil.
7. Carefully peel the electrodes from the clear plastic film - Do not throw this plastic film away.
8. Attach the electrode(s) to the appropriate muscle and insert the PERIFORM®+ or ANUFORM® if you are using it. The reference electrode should be placed over a bony point i.e. for pelvic floor retraining place on the hip area.
9. Turn the unit on by pressing ON/OFF button once.

10. Set up the Threshold level microvolt (μV). This is achieved by contracting your muscles for 5 seconds and taking a reading of the peak value. The values on the LCD screen will be moving up and down. Take approximately 40% of the peak value. E.g. if the reading moves from 25-30 μV , take 40% of 30 μV which is 12 μV as the threshold level. This level will be different for each individual.
11. Set the Threshold level by pressing the + or - button.
12. Set the desired parameters by pressing the SET button. This will move you through the parameters for each mode (VOL, WRK, RST, TRS, FAB, DOF, WDE – see page 12 for further details). Press the + or - buttons to change each parameter.
13. Place the unit on a desk in its stand or hold the unit.
14. Relax so that the microvolts reading is as low as you can manage.
15. To start the program, press the START/ESC button. To end the program before it is complete, press the START/ESC button again or turn the unit off.
16. When you have finished, remove and replace the electrodes onto the clear plastic film and reseal them in the plastic bag. If you are using a PERIFORM[®]+ or ANUFORM[®], clean the probe and let it dry before sealing it in its plastic bag. See cleaning information on page 19.
17. Store the surface electrodes in a cool place, such as a fridge to extend their 'in-use' life.

NOTE: If the parameters of the unit cannot be changed when pressing the SET button, the unit will need to be unlocked. Use the pin from one of the lead wires and press it into the concealed LOCK button at the back of the unit

Advanced Instructions

Firstly, the EMG Threshold level (THRS) requires setting up as accurately as possible. The threshold level will vary from one patient to another. The threshold level is a goal or target a patient strives to achieve when contracting their muscles or when performing exercises.

To set up the threshold the patient should be asked to contract and/or perform the exercise holding the position for approx. 5 seconds and then relax for 5-10 seconds before repeating the contraction or movement. Note the peak microvolt μV reading displayed on the LCD. The microvolt readings will change due to muscle variation.

Calculate the average of the two peak readings and set the threshold between 40-50% of this reading.

Example: Average Peak reading: $30\mu\text{V}$ take 40% of this value = $12\mu\text{V}$ threshold setting.

To adjust the threshold setting (shown at the top of the LCD screen) press the THRS + or – buttons.

Setting the EMG parameters

- Pressing the SET button consecutively takes you through a looped menu from which you can change the parameters for biofeedback including volume, work/rest times, trials, feedback sound, serial data and band settings.

Press SET to see:	
VOL	Press the THRS + or THRS – button to adjust the sound volume between 0 – 10.
WRK	Press the THRS + or THRS – button to adjust the work time in seconds (2-99).
RST	Press the THRS + or THRS – button to adjust the rest time in seconds (2-99).
TRS	Press the THRS + or THRS – button to adjust the number of trials (repetitions).
FAB FBL FOF	Press the THRS + or THRS – button to select either: FAB - to get audio feedback when the patient exceeds the threshold. FBL - to get audio feedback when patient goes below the threshold. FOF - to turn off any audio biofeedback.
DON DOF	Press the THRS + or THRS – button to switch serial data on or off. Switch on if connected to PC (this will stop the device from turning off).
WDE NRW	Press the THRS + or THRS – button to select WDE wide filter band or NRW narrow filter band.

Electrodes placed on the legs, arms, face buttock or pelvic area use WIDE BAND setting. All other areas i.e. abdomen, chest, shoulder, upper arms and back use the NARROW BAND setting.

Narrow band setting is used due to the interference from the electrical activity of the heart which needs to be eliminated from the EMG signal.

Pressing the ESC button will record and store the parameters that have just been set until it is reprogramed. To start the program press the START button, to escape press the ESC button.

At the end of each session for example Work 5 seconds, Rest 5 seconds and 5 trials (Repetitions) the LCD display will highlight the statistics (described on page 15 and 16). The results will differ from one patient to another.

EMG Statistics.

Once a session has finished WAV will appear on the screen. Press the SET button consecutively to read all the statistics.

Press SET to see:	
WAV	This is the work average for the session measured in (μ V) microvolts. Average readings will vary from one patient to another.
RAV	This is the rest average for the session measured in μ V - Below 4 μ V a muscle is beginning to rest.
OST	This is the average onset of muscle contraction measured in seconds, readings below 1 second can be considered normal.
RLS	This is the average muscle release measured in seconds, readings below 1 second can be considered normal.

When you have completed reading the statistics you can either continue to press the SET button to read the statistics again or press the ESC button to come out of the program. NOTE: Once you have pressed the ESC button you will lose all the information.

If you are linked to a PC the statistics can be saved for future reference.

Recording Information

To record and store the information for one patient over a number of sessions using the LOCK function mode


- Remove the battery compartment cover.
- Using a thin rod (end of the lead wires) press the concealed LOCK button once (see diagram on page 9) and you will hear two beeps, this indicates the unit is locked. A padlock icon will appear on screen.
- A series of sessions can then be undertaken. After each session press the SET button consecutively to read the statistics for that particular session. After reading all the statistics, press the ESC button. The information at this point will be lost for good.
- To view the statistics over a number of sessions press the concealed LOCK button once to unlock the unit. You will hear one beep. All session data will have been added together to form the overall average reading for the given period. This includes the number of trials performed and average muscles activity.
- Press the SET button consecutively to view these statistics.
- After reading all the statistics, press the ESC button. The information at this point will be lost for good.
- Don't forget to LOCK the unit again if you wish to record the average statistics of another series of sessions.

If you are linked to a PC read the software instructions for downloading patient recorded details.

NOTE: Programs that are not completed will not be recorded.

NOTE: While the unit is in operation only the volume and threshold can be changed. Recording only occurs by pressing the start button and when the program is completed.

Statistics


EXAMPLE

The Work / Rest statistics automatically shown at the end of every Work / Rest cycle are derived as follows:

1. **TLS** - Trials: Displays the number of work rest trials for the entire period the unit was locked.
2. **WAV** - Work Average: The average value in microvolts of all the Work segments excluding the first second of each segment.
3. **RAV** - Rest Average: The average value of the Rest segments excluding the first second of each segment.
4. **OST** - Average Onset: The average time taken after each respective "Work" prompt for the EMG to reach 75% of the average value of the third second of that particular work segment. In the case that only two seconds of work has been selected then, Average Onset is the average of the time taken to reach 75% of the average value of the second second (after the "Work" prompt). If any of the Onsets in any of the trials is more than 2 seconds then it (or they) will be rejected and indicated by a dash before the number displayed. The display will indicate only the average of those trials which were 2 seconds or less. If all the Onsets were more than 2 seconds then the display will read AB2.

(continued)

5. **RLS - Average Release:** The average time taken, after the REST prompt to reach 37.5% of the previous third second of work or as in item 4 above, the second second if only 2 seconds of work has been selected. As in the case of Onset, if any Release of the trials is more than 2 seconds then it (or they) will be rejected and the display will indicate the average of only those which were 2 seconds or less but a dash (-) before the number displayed will indicate this rejection.
6. If all the Releases were more than 2 seconds then the display will read AB2.

Electrode Types & Tips

Both internal and external electrodes can be used with this unit.
They are single patient use only.

Internal electrodes:

021972 PERIFORM®+ (vaginal electrode) 011503 ANUFORM® (anal electrode)

External self adhesive electrodes:

012571 50x50mm

012572 90x50mm

012574 30mm diameter

012575 50mm diameter

Tips

- If using an internal electrode, after use always wash it in warm soapy water, rinse and dry thoroughly and carefully store it for your next treatment session.
- Never let anyone else use your electrode.
- Follow all instructions provided with the electrode.
- If using self adhesive electrodes always cleanse the skin with soap and water, then rinse and dry well in order to ensure good adhesion of the electrodes.
- Clip away hairy skin using scissors. Do not shave the area.
- After using surface electrodes, place them back on the plastic film they were packaged in and store them in a cool place i.e. the fridge.
- After repeated use, the sticky surface on the electrode will eventually become dry. This surface is water based so by adding a few drops of water will give you a few extra days of electrode life back!
- If you find the electrodes will not stick due to oily skin, cleanse the skin with soap and water, then rinse and dry the area around the electrode site. If this does not work, try cleansing the skin with a swab impregnated with alcohol.
- The electrodes conductive material is water- based. If it becomes saturated (e.g. from perspiration), it will lose its adhesive qualities. After use leave the electrodes face up overnight to dry out (replace on plastic film in the morning). At some point the electrodes will become dry. Moisten the adhesive surface with a few drops of water, and apply onto the plastic film overnight. This procedure will increase the electrode life by few more days.

Care and Maintenance

Control Unit:

- Wipe the surface after use with a damp cloth or antiseptic wipe.
- Do not use cleaning sprays or alcohol based cleaning solutions.

Battery:

- To change the battery, open the battery door on the rear of the control unit by pressing down on the raised rib pattern just below the belt clip. Lift the battery out of the compartment. This is very easy and can be done by the user.
- Check periodically for any discharge from the battery.
- Remove battery completely from unit if not in use for any extended period of time (typically one week).
- Low battery indicator of 6.9 volts shown on LCD display. When flashing change battery for a new one.
- Preferably use a PP3 alkaline battery.

Lead Wires:

- The lead wires should be handled carefully and never stretched, as this can cause the stimulation to function below normal standards or not at all.
- Examine lead wires before each treatment for loose connections or damage.
- Avoid twisting the lead wires.
- Wipe lead wires after use with a damp cloth or antiseptic wipe.
- Store the lead wires carefully after each use.

Self-Adhesive Electrodes:

- Check the short connectors have not become separated from the electrodes.
- Replace electrodes onto plastic film after use. If they drop onto the floor debris will adhere to conductive gel rendering the electrodes ineffective.

Vaginal / Anal Electrodes:

- Check the connectors have not become separated from the probe.
- Before & after use always wash it in warm soapy water, rinse and dry and carefully store.
- Never let anyone else use your electrode.
- Follow all instructions provided with the electrode.

Caution: Static electricity may damage this product.

NOTE: Only Patterson Medical or appointed distributors/importers are approved to undertake servicing.

Specifications

1. Single Channel
2. EMG Range: 0.2 to 2000 μ V RMS (continuous)
3. Sensitivity: 0.1 μ V RMS
4. Accuracy: 4% of μ V reading \pm 0.3 μ V at 200 Hz
5. Selectable Bandpass filter - 3db Bandwidth;
 - a. Wide: 18 Hz \pm 4 Hz to 370 Hz \pm 10% - Reading below 235 microvolts. 10 Hz \pm 3 Hz to 370 Hz \pm 10% - Reading above 235 microvolts.
 - b. Narrow: 100 Hz \pm 5% to 370 Hz \pm 10%
6. Notch filter: 50 Hz (Canada 60Hz) - 33 dbs (0.1% accuracy)
7. Common Mode Rejection Ratio: 130 dbs Minimum @ 50 Hz
8. Battery: PP3 Alkaline
9. Work / Rest periods: 2-99 seconds
10. Number of Trials: 1-99
11. Low battery indication at 7.4 volts \pm 0.2 volts and auto shut off 10 minutes after last key pressed
12. Environmental conditions for storage & transport:
-10 to +50°C. 0-90% Humidity
13. Physical Dimensions: Length 128.5mm, Width 64mm, Depth 28.3 mm
14. Weight: 0.15 Kg

Interference

The PERITONE may be subject to Electromagnetic Interference (see page 5).

Additionally, the power supplies of some notebook computers can give off substantial amounts of interference to which the PERITONE is susceptible. This is particularly the case when the power supply “block” is the type with only a two pin connector connecting it to the mains with no earth.

As a precaution, make sure that the power cable from the notebook runs as far away as possible from the patient connection wires of the PERITONE. (This does not include the fibre optic cable as it is, electrically, entirely non-conductive).

Try to keep the PERITONE close to the patient’s body (in the “field” of the patient) either on their lap, in their pocket or clipped to their belt. Also run the electrode wires along the surface of the patient’s body keeping the wires as close as is practically possible all the way from the electrodes to the PERITONE unit and not dangling freely in space away from the patient.

A relaxed muscle should read below 3.5 microvolts (μV). If you find that, even when the patient’s muscle is soft and relaxed to the touch, you are still reading persistent, unexpectedly high μV levels, try turning off the notebook’s external mains power supply. (The notebook will continue to run on it’s own internal battery). If the μV reading(s) suddenly reduce and then go back up when you turn on the notebook’s power supply then this is the cause of the interference. If this problem causes ongoing difficulties, please contact Patterson Medical for assistance.

Troubleshooting

If you experience problems with EMG readings we recommend you make the following checks:

1. Check the lead wires for splits or breaks in the wire or at the end where the connectors are attached to the wire.
2. Check the lead wires of surface and or internal electrodes. Inferior surface electrodes will cause incorrect readings; we recommend you always use quality electrodes for EMG measurement.
3. If you are using a laptop computer and experience interference when using the charger, switch the charger off. If you still experience interference move to another area.
4. If you are using vaginal or anal probes we recommend the patient uses conductive gel.
5. Some patient's vaginal aperture may be too large for some internal probes, causing intermittent contact with the walls of the pelvic muscle. In such cases, one should try another larger electrode.
6. If you are connected to a laptop or desk top computer check to make sure the ends of the fibre optic cable are not damaged, as this will tend to obstruct the signal from the PERITONE unit to the computer.
7. If you are not connected to a computer check to make sure there is no debris in the fibre optic connector, which is situated on the top left hand side of the device.

Software

An optional computer program is available for use with the PERITONE unit. It connects to the computer via a fibre optic cable and serial dongle.

The software can be purchased separately from Patterson Medical.

The software displays the EMG results in the form of a graph and has the following modes associated with it:

- Open display
- Template training
- Work / Rest assessment
- Work / Rest training

Single channel progress and test reports can be created from the work / rest assessments. These enable the user to view up to 10 assessments on the computers screen or in a printed report over a period of time. These reports can then be used for insurance purposes or to prove that the patient is improving.

Please contact your distributor or Patterson Medical for more information.

Software

File Reports **Graph** Setup Help

- Print patient details
- Print patient sessions
- Print patient list
- Print patient session statistics
- Create patient progress report
- Saved progress reports
- Load the EMG Real Data CSV file

Select Patient

Enter Patient ID or Surname to search: _____ Session Condition: _____ Incontinence

Progress Report Channel A Only
Test Report Channel A and Channel B

Patient	Surname	First name
1	Lennon	John
1A	Example	Bob
2	Purple	Andy
3	Purple	Freddie

☐ Confidential ☐ Default Patient


Patient Options: View / Edit View Daily Report Add Download Delete Patient Record


Session Options: Review Notes Delete Session ☐ Click on a checkbox to Select or De-Select sessions for a Progress report. De-select All

Close


S	Date	Time	Protocol	Condition
<input type="checkbox"/> 2	20/03/2008	17:15	Open Display	Incontinence
<input type="checkbox"/> 3	20/03/2008	17:17	ETS	Incontinence
<input type="checkbox"/> 4	20/03/2008	17:18	Templates	Incontinence
<input checked="" type="checkbox"/> 5	20/03/2008	17:19	Assessment 51	Incontinence
<input checked="" type="checkbox"/> 6	20/03/2008	17:20	Assessment 51	Incontinence
<input checked="" type="checkbox"/> 7	20/03/2008	17:22	Assessment 51	Incontinence

PC Software: Patient Database PC Software: Clinical Progress Report

PC Software: Template training - draw the Template shape, multiply it and follow by your EMG.

Disposal

Information on disposal for users of Waste Electrical & Electronic Equipment (WEEE) for private households:

Electrical and electronic products including batteries should not be mixed with general household waste. For proper treatment, recovery and recycling, please take these products to designated collection points, where they will be accepted on a free basis.

Applications

- Promotes continence
- Increases muscle strength
- Maintains or improves range of movement
- Increases and improves the blood supply to the muscle
- Reduces pain

Information regarding Electromagnetic compatibility and interference (EMC)

PERITONE is designed to produce very low levels of radio frequency (RF) emissions (interference), to be immune from effects of interference produced by other equipment operating in their vicinity and damage due to electrostatic discharge all when operating in a typical domestic and or clinical environment. PERITONE is certified to meet the international EMC standard EN60601-1-2. For more information please refer to the tables 201, 202, 204 and 206.

In the EMG mode the PERITONE may be subjected to Electromagnetic Interference.

Additionally, the power supplies of some notebook computers can give off substantial amounts of interference which the PERITONE is susceptible.

This can happen when the power supply "block" has only a two pin connector connecting it to the mains with no earth.

As a precaution, make sure that the power cable from the notebook is placed as far away as possible from the connection wires of the PERITONE. (This does not refer to the fibre optic cable which is electrically non-conductive).

Try to keep the PERITONE close to the patient's body (in the "field" of the patient) either on their lap, in their pocket or clipped to their belt. Keep the electrode wires as close as possible to the patients' body and not dangling freely.

A relaxed muscle should read below 3.5 microvolts (μV). If even when the patient's muscle is soft and relaxed to the touch, the reading is still high, try turning off the notebooks external main power supply. (The notebook will continue to run on its own internal battery). If the μV reading(s) suddenly reduce(s) and then go back up after turning on the notebook power supply, it means that an interference has occurred.

**Table 201: Guidance and manufacturer's declaration
– electromagnetic emissions**

The PERITONE product is intended for use in the electromagnetic environment specified below. The customer or the user of the PERITONE product should ensure that it is used in such an environment.		
Emissions test	Compliance	Electromagnetic environment – guidance
RF emissions CISPR 11	Group 1	The PERITONE product uses RF energy only for its internal function. Therefore, its RF emissions are very low and are not likely to cause any interference in nearby electronic equipment.
RF emissions CISPR 11	Class B	The PERITONE product is suitable for use in all establishments, including domestic establishments and those directly connected to the public low-voltage power supply network that supplies buildings used for domestic purposes.
Harmonic emissions IEC 61000-3-2 IEC 61000-3-2	Not applicable	
Voltage fluctuations/ flicker emissions IEC 61000-3-3	Not applicable	

**Table 202: Guidance and manufacturers declaration
– electromagnetic immunity**

The PERITONE product is intended for use in the electromagnetic environment specified below. The customer or the user of the PERITONE product should ensure that it is used in such an environment, and that precautions regarding that environment are heeded.			
Immunity test	IEC 60601 test level	Compliance level	Electromagnetic environment – guidance
Electrostatic discharge (ESD) IEC 61000-4-2	±6 kV contact ±8 kV air	±6 kV contact ±8 kV air	Floors should be wood, concrete or ceramic tile. If floors are covered with synthetic material, the relative humidity should be at least 30 %.
Power frequency (50/60 Hz) magnetic field IEC 61000-4-8	3 A/m	3 A/m	Power frequency magnetic fields should be at characteristic levels of a typical location in a typical commercial or hospital environment.

Table 204: Guidance and manufacturer's declaration – electromagnetic immunity

<p>The PERITONE product is intended for use in the electromagnetic environment specified below. The customer or the user of the PERITONE product should ensure that it is used in such an environment.</p>			
Immunity test	IEC 60601 test level	Compliance level	Electromagnetic environment – guidance
<p>Conducted RF</p> <p>IEC 61000-4-6</p> <p>Radiated RF</p> <p>IEC 61000-4-3</p>	<p>3 Vrms</p> <p>150 kHz to 80 MHz</p> <p>3 V/m</p> <p>80 MHz to 2,5 GHz</p>	<p>3 Vrms</p> <p>150 kHz to 80 MHz</p> <p>3 V/m</p> <p>80 MHz to 2,5 GHz</p>	<p>Portable and mobile RF communications equipment should be used no closer to any part of the PERITONE product, including cables, than the recommended separation distance calculated from the equation applicable to the frequency of the transmitter.</p> <p>Recommended separation distance</p> $d = 1.2\sqrt{P} \text{ (150 kHz to 80 MHz)}$ $d = 1.2\sqrt{P} \text{ (80 MHz to 800 MHz)}$ $d = 2.3\sqrt{P} \text{ (800 MHz to 2.5 GHz)}$ <p>where P is the maximum output power rating of the transmitter in watts (W) according to the transmitter manufacturer and d is the recommended separation distance in meters (m).</p> <p>Field strengths from fixed RF transmitters, as determined by an electromagnetic site survey,</p> <p>(a) should be less than the compliance level in each frequency range;</p> <p>(b) interference may occur in the vicinity of equipment marked with the following symbol: </p>
<p>NOTE 1: At 80 MHz and 800 MHz, the higher frequency range applies.</p> <p>NOTE 2: These guidelines may not apply in all situations. Electromagnetic propagation is affected by absorption and reflection from structures, objects and people.</p>			
<p>(a) Field strengths from fixed transmitters, such as base stations for radio (cellular/cordless) telephones and land mobile radios, amateur radio, AM and FM radio broadcast and TV broadcast cannot be predicted theoretically with accuracy. To assess the electromagnetic environment due to fixed RF transmitters, an electromagnetic site survey should be considered. If the measured field strength in the location in which PERITONE product is used exceeds the applicable RF compliance level above, the PERITONE product should be observed to verify normal operation. If abnormal performance is observed, additional measures may be necessary, such as reorienting or relocating the PERITONE product.</p> <p>(b) Over the frequency range 150 kHz to 80 MHz, field strengths should be less than 3 V/m.</p>			

Table 206: Recommended separation distances between portable and mobile RF communications equipment and PERITONE product

The PERITONE product is intended for use in an electromagnetic environment in which radiated RF disturbances are controlled. The customer or the user of the PERITONE product can help prevent electromagnetic interference by maintaining a minimum distance between portable and mobile RF communications equipment (transmitters) and the PERITONE product as recommended below, according to the maximum output power of the communications equipment.

Rated maximum output power of transmitter W	Separation distance according to frequency of transmitter		
	150 kHz to 80 MHz	80 MHz to 800 MHz	800 MHz to 2,5 GHz
	$d = 1.2\sqrt{P}$	$d = 1.2\sqrt{P}$	$d = 2.3\sqrt{P}$
0,01	0.12	0.12	0.23
0,1	0.38	0.38	0.73
1	1.2	1.2	2.3
10	3.8	3.8	7.3
100	12	12	23

For transmitters rated at a maximum output power not listed above, the recommended separation distance d in meters [m] can be estimated using the equation applicable to the frequency of the transmitter, where P is the maximum output power rating of the transmitter in watts (W) according to the transmitter manufacturer.

NOTE 1: At 80 MHz and 800 MHz, the separation distance for the higher frequency range applies.

NOTE 2: These guidelines may not apply in all situations. Electromagnetic propagation is affected by absorption and reflection from structures, objects and people.

Warranty

Patterson Medical provides a warranty to the original purchaser that this product will be free from defects in the material, components and workmanship for a period of 2 years from the date of purchase.

If Patterson Medical is satisfied that the product/s is defective the purchaser may return this unit/s to Patterson Medical or the appointed distributor for repair or replacement with a new unit.

All returns must first be authorised by Patterson Medical in advance. The liability of Patterson Medical under this limited product warranty does not extend to any misuse or abuse such as dropping or immersing in water or other liquid substance or tampering with the unit or normal wear and tear. Any evidence of tampering will nullify this warranty.

Customer Service

If you have any questions contact our Customer Service Team:

Patterson Medical Ltd. Nunn Brook Road, Huthwaite,

Sutton-in-Ashfield, Nottinghamshire, NG17 2HU. UK

Tel: +44 (0)8448 730 035 Fax: +44 (0)8448 730 100

Email: mobilisrolyan.sales@patterson-medical.com

Web: www.mobilisrolyan.com

The Neen PERITONE is CE Approved

CE0088

Manufactured under a quality management system in accordance with ISO9001 and ISO13485 by LRQA 0964891

Clinical References

Wilke, H.J., Wolf, S., Claes, L.E et al (1990); Stability increase of lumbar spine with muscle groups. A biomedical in vitro; Spine 20,192-198. P31

Roy, S.H., De Lueas, C.J., Casavant, D.A. (1992); Shoulder torque changes in the swimming athlete; Am. J. Sports Med. Rehab. 20, 323-335.

Lucca, J.A., Recchiutu, S.J. (1983); Effect of electromyographic biofeedback on an isometric strengthening program; Phys. Therapy 63, 200-203.

Nieuwmeyer, P.A., Visser, S.L., Feenstral, L. (1985); Bell's palsy: A Polyneuropathy; Am. J. Otol. 6, 250-252.

Booker, H.E., Rubow, R.T., Coleman, P.J. (1996); Simplified feedback in neuromuscular retraining. An automated approach using electromyographic signals; Arch. Phys. Med. Rehab. 50, 621-625.

Reld, D.C., Saboe, L.A., Chepan, J.C. (1996); Shoulder instability, biofeedback, isokinetics, evidence-based practice; Physiotherapy Canada 1, 251-256.

Cardozo, L.D., Stanton, S.L., Hafner, J., Allan, V. (1978); Biofeedback in the treatment of detrusor instability; British Journal of Urology 50, 250-254.

Sherman, R., Arena, J., Sherman, C., Ernst, J. (1989); The mystery of phantom pain: Growing evidence for physiological mechanisms; Biofeedback and self-regulations 14, 267-280.

Gowen, I., Jobe, F., Tibone, J. (1987); A comparative electromyographic analysis of the shoulder during pitching; Am. J. Sports Med 50, 586-590.

McConnell, J.S. (1986); The management of chondromalacia patella: A long term solution; Australian J. Physiotherapy 32 [4], 215-223.

Doham, R.S., Lakin, D.M. (1978); An evaluation of electromyographic biofeedback in the treatment of myofacial pain-dysfunction syndrome; J. Am. Dent. Assoc., 96: 656-662.

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Logiciel

File Reports Graph Setup Help

Print patient details
Print patient sessions
Print patient list
Print patient session statistics
Create patient progress report
Saved progress reports
Load the EMG Real Data CSV file

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☐ Confidential ☐ Default Patient


Patient Options



View / Edit View Daily Report
Add Download
Delete Patient Record

Session Options

Review Click on a checkbox to Select or De-Select sessions for a Progress report
Notes
Delete Session De-select All

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